

MICROCT DATA PROCESSING GUIDELINES

Andy Kiss

August 24, 2017

ABSTRACT

Here is a brief guideline for processing microtomography data collected using the microCT at SLAC/SSRL. This document can be used to reconstruct projection data into meaningful 3D datasets but it should be noted that some tweaks may need to be performed from this outline.

CONTENTS

Abstract.....	2
Contents.....	2
Prerequisites	3
Convert TIFF Stacks to Reference Corrected Projection Images	4
Mosaic Stitching.....	6
Find the Rotation Center.....	7
Create Sinograms	8
Sinogram Reconstruction.....	9
Opening the Reconstructed Slices	10
Python	10
Fiji.....	10
TXM Wizard.....	11


PREREQUISITES

Below is a list of software used in the following outline.


- Python 3+
 - Anaconda distribution (recommended), using Spyder
 - Go to website <https://www.continuum.io/downloads> and download Anaconda installation package.
 - Packages:
 - Scipy
 - Numpy
 - Pywavelets
 - ASTRA (<http://www.astra-toolbox.com>)
- Py_image_processing
 - A series of scripts on GitHub for this reconstruction
 - These can be downloaded from https://github.com/slaclab/py_image_processing
- FIJI (ImageJ 2.0)
 - Not necessary but helpful for viewing and manipulating data
 - Go to website <https://fiji.sc/#download> and download FIJI zipped folder.
 - Unzip the FIJI folder
 - In Fiji.app folder double click ImageJ-win32 to launch the software.
- TXM Wizard
 - Necessary for mosaic imaging
 -

CONVERT TIFF STACKS TO REFERENCE CORRECTED PROJECTION IMAGES

- Locate the folder where you saved your data. There should be a scan log file (*.log) as well as folders with your transmission images, dark field, and flat field image stacks.



- Open Spyder and open the *make_projections.py* script.
- Locate the lines that set the scan log file name, *fn_scanlog*, and path, *path*, and set the variables to point to your data.



- Scan through the remaining settings in that block. Many of them can stay the same from reconstruction to reconstruction. If you are interested in modifying the binning of the images or removing zingers, here is where the settings are.
- Note: if you have a data processing computer with sufficient memory to perform all the analysis in memory (>32 GB), check that *flag_big_memory* is set to *True*. This will run the calculations much faster.
- Note: By default, the program will export *.bim files. These are binary image files with the metadata stored in the header. It is recommended that you perform the analysis using this file format. Check that *flag_bim* is set to *True*, otherwise it will export *.tif files.
- When ready, save and run make_projections.py. This will take a while but there will be output to the console so you can track the process.
- The script will make a post-processing folder, *sample-name-processing*, where it will store the averaged dark and flat field images, and a folder where it will export the reference corrected projections.

The screenshot shows the Spyder Python IDE with the file `make_projections.py` open in the editor. The script contains various parameters for image processing, including file paths, binning settings, and flags for averaging and memory usage. The IPython console at the bottom shows the execution progress, with messages like "Working on image (000708/000721)...done" repeated for each image in the dataset. The console also shows the final output "Done." and the prompt "In [2]:".

```

1 # -*- coding: utf-8 -*-
2 """
3 Created on Wed Jan 25 22:39:41 2017
4
5 @author: Andy
6 """
7
8
9 import os
10 import numpy as np
11 import tsm_image
12 import image_handling
13
14
15 # %% User needs to set these parameters
16 fn_scanlog = '20170330_153537_LLNL_BadFi.log'
17 path = r'C:\Users\andykiss\Documents\recon_demo'
18
19
20 # Despeckling images
21 flag_despeckle = True
22 DS_delta = 500
23 DS_rad = 3
24
25 # Bin images
26 flag_bin = False
27 B = 2 # Amount to bin (ex. B = 2 --> 2x2 binning)
28
29 # Exposure times
30 t_ref = 500 # [ms] reference exposure time
31 t_proj = 500 # [ms] projection exposure time
32
33 # True -- Average any images together
34 # False -- Use median averaging
35 flag_average = True
36
37 # Left-right mosaic so that the center of sample rotation is at the edge of the
38 # field of view so 180 deg is equivalent to 0 flipped and translated
39 flag_360_mosaic = False
40 new_xy = (500., 0.)
41
42 # True -- Load everything into memory, and then write to file. Faster.
43 # False -- Load files and write one by one. Slower.
44 flag_big_memory = True
45
46 # True -- Flatten one of the big mosaic images (image stack 000708-000721)

```

Name	Type	Size	Value
B	int	1	2
DS_delta	int	1	500
DS_rad	int	1	3
N	int	1	2
N_proj	int	1	721
Nh	int	1	360
e	float	1	2.718281828459045
euler_gamma	float	1	0.5772156649015329
flag_360_mosaic	bool	1	False

IPython console

```

Working on image (000708/000721)...done
Working on image (000709/000721)...done
Working on image (000710/000721)...done
Working on image (000711/000721)...done
Working on image (000712/000721)...done
Working on image (000713/000721)...done
Working on image (000714/000721)...done
Working on image (000715/000721)...done
Working on image (000716/000721)...done
Working on image (000717/000721)...done
Working on image (000718/000721)...done
Working on image (000719/000721)...done
Working on image (000720/000721)...done
Working on image (000721/000721)...done
Done.
In [2]:

```

Permissions: RW End-of-lines: CRLF Encoding: UTF-8 Line: 34 Column: 32 Memory: 8 %

MOSAIC STITCHING

TXM Wizard

- Open *Mosaic Stitcher*
- Click *Browse* and load only the images necessary for the first mosaic image
 - All of the imaging properties should load into the fields on the right
- Start to use the options explained below to stitch together those. Click *Run Stitching* to see the stitched image.
 - *Use PhCo Alignment* will use phase correlation to help align the images for stage errors
 - *Normalization among Tiles* will normalize the edges
 - *Smooth The Edge* will blend the images together
- Once satisfied with the result, enable *MultipleCPU* and make sure the images are saved as **.bim*
- Click *Batch Stitching* and select all the images

FIND THE ROTATION CENTER

- In Spyder, open *rotation_center.py*
- Find the Registration Settings cell and change
 - *root* to point to the folder with the reference corrected *.bim files
 - *fn0* to the file name of the first image in your dataset (usually at theta = 0°)
 - *fn1* to the file name of the image in the dataset rotated 180° from the first file
- Save and run the script. The console will output the X and Y shift, as well as the pixel location of the rotation axis.

The screenshot shows the Spyder Python IDE with the *rotation_center.py* script open in the editor. The script includes comments and code for finding the rotation center. The Variable explorer on the right shows the state of variables after execution. The IPython console at the bottom displays the output of the script, including the registration settings and the calculated offset and rotation axis.

Name	Type	Size	Value
B	int	1	2
D5_delta	int	1	500
D5_rad	int	1	3
I0	float32 (2048, 2048)		array([[0.05339661, 0.06957942, 0.0280 ...
I1	float32 (2048, 2048)		array([[0.02074633, 0.02706698, 0.0 ...
I1_reg	float32 (2048, 2048)		array([[0.02660256, 0.02018541, 0.0 ...
I_diff	float32 (2048, 2048)		array([[0.02679406, 0.04939402, 0.0280 ...
N	int	1	2
N_proj	int	1	721


IPython console output:

```
In [2]: runfile('C:/Users/andykiss/Documents/programming/python/py_image_processing/rotation_center.py', wdir='C:/Users/andykiss/Documents/programming/python/py_image_processing')
Reloaded modules: image_handling.resize_image, image_handling.wavelet_fourier_filter, tsm_image.microCT_scanlog, image_handling.average_image, tsm_image.formats.binsino, image_handling, tsm_image, image_handling.remove_outliers, tsm_image.add_metadata, tsm_image.formats, tsm_image.formats.binfile, tsm_image.formats.bim, image_handling.generate_synthetic_data, image_handling.formats.tiff

The offset is (deltx, dely) = (6.10, 0.20)
The rotation axis is at 1017.90

In [3]:
```

- A figure will show the loaded images and the difference map after they have been aligned. The quality of the alignment can be judged based on the difference map. If it is easy to see the features in the original images, then the alignment is poor and will not lead to a high quality reconstruction.



CREATE SINOGRAMS

- In Spyder, open `create_binsinos.py`
- Find the Settings cell and change the path variable, `root`, to the folder with the reference corrected projections
- Set the rotation axis, `rot_axis`, to the rotation axis value found in the previous step
- Ring artifacts in the reconstruction can be removed by using a filter. Set the ring removal filter flag, `flag_ring_removal`, to `True` to remove these.
- Run the script and wait for the sinograms to be output to the `sinos` folder

The screenshot shows the Spyder Python IDE interface. The main editor displays the `create_binsinos.py` script. A red box highlights the settings section of the script, which includes:

```
22 # %% Settings
23 # File names and Location
24 root = r'C:\Users\andykiss\Documents\recon_demo\20170330_153537_LLNL_BadTi-proc
25 ext = '.bim'
26
27 # True rotation axis
28 rot_axis = 1017.9
29
30 # Ring removal settings
31 flag_ring_removal = True
32 wf_str = 'sym16'
33 wf_n = 6
34 wf_sig = 2.0
35
36 # Output directory
37 outdir = 'sinos/'
38
39
40 # %% Find the files
41 if (root[-1] != '\\' and root[-1] != '/'):
42     root += '/'
43
44 ls = os.listdir(root)
45 for fn in ls:
```


The 'Variable explorer' panel on the right shows the following variables:

Name	Type	Size	Value
path_process	str	1	20170330_153537_LLNL_BadTi-processing/
path_proj	str	1	20170330_153537_LLNL_BadTi_1/
path_ref1	str	1	20170330_153537_LLNL_BadTi-ref_1/
path_ref2	str	1	20170330_153537_LLNL_BadTi-ref_2/
pi	float	1	3.141592653589793
root	str	1	C:\Users\andykiss\Documents\recon_d...
rot_axis	float	1	1017.9
row	int	1	2048
shift	int	1	12

The IPython console at the bottom shows the output of the script, indicating that sinograms are being written for various years (2034/2048) and that the process is complete.

SINOGRAM RECONSTRUCTION

- In Spyder, open *ASTRA_recon_batch.py*
- Set the path to the sinograms, *path*
- Define your reconstruction settings
 - If possible, use a CUDA enabled algorithm as the GPU is much faster than the CPU
- Run the script



The screenshot shows the Spyder Python IDE interface. The main editor displays the `ASTRA_recon_batch.py` script. A red box highlights the section where the path and reconstruction settings are defined. The IPython console at the bottom shows the execution progress, listing files being reconstructed from 0005/2048 to 0024/2048.

```
# -*- coding: utf-8 -*-
"""
Reconstruct the *.binsino files in a folder using ASTRA
Created on Fri Aug 18 10:04:37 2017
@author: andykiss
"""
import os
import tom_image
from tomo_recon import astra_recon

# %% Define the folders
path = r'C:\Users\andykiss\Documents\recon_demo\20170330_153537_LLNL_BadTi-process'
ext = '.binsino'

# ASTRA settings
alg = 'FBP_CUDA'
alg_iter = 100
px = 1.0
flag_matlab_crop = True

# %% Get the files
if (path[-1] != '\\') and path[-1] != '/':
    path += '/'
ls = os.listdir(path)
for fn in ls:
    if (not fn.endswith(ext)):
        ls.remove(fn)
N = len(ls)

# %% Start reconstructing
# Make output directories
outdir = 'slices/'
if (not os.path.isdir(path + outdir)):
    os.mkdir(path + outdir)
outdir = outdir + alg + '/'
if (not os.path.isdir(path + outdir)):
    os.mkdir(path + outdir)
```

Name	Type	Size	Value
mod_path	str	1	C:\Users\andykiss\Documents\program...
new_xy	tuple	2	(500.0, 0.0)
outdir	str	1	slices/SIRT_CUDA/
path	str	1	C:\Users\andykiss\Documents\recon_d...
path_df	str	1	20170330_153537_LLNL_BadTi-df_1/
path_process	str	1	20170330_153537_LLNL_BadTi-pr...
path_proj	str	1	20170330_153537_LLNL_BadTi_1/
path_ref1	str	1	20170330_153537_LLNL_BadTi-ref_1/
path_ref2	str	1	20170330_153537_LLNL_BadTi-ref_2/

IPython console

```
Reconstructing file (0005/2048)...done
Reconstructing file (0006/2048)...done
Reconstructing file (0007/2048)...done
Reconstructing file (0008/2048)...done
Reconstructing file (0009/2048)...done
Reconstructing file (0010/2048)...done
Reconstructing file (0011/2048)...done
Reconstructing file (0012/2048)...done
Reconstructing file (0013/2048)...done
Reconstructing file (0014/2048)...done
Reconstructing file (0015/2048)...done
Reconstructing file (0016/2048)...done
Reconstructing file (0017/2048)...done
Reconstructing file (0018/2048)...done
Reconstructing file (0019/2048)...done
Reconstructing file (0020/2048)...done
Reconstructing file (0021/2048)...done
Reconstructing file (0022/2048)...done
Reconstructing file (0023/2048)...done
Reconstructing file (0024/2048)...done
```


Permissions: RW End-of-lines: CRLF Encoding: UTF-8 Line: 24 Column: 11 Memory: 10 %

OPENING THE RECONSTRUCTED SLICES

- The reconstructed slices can be opened using Python, Fiji, or TXM Wizard



Python

- Import the txm_image package
- Use the read_file function to load the image



Fiji

- File -> Import -> Raw...
- Choose the file you would like to open using the file dialog. If you want to open all the files in the folder, still only choose one.
- Use the settings listed below, 32-bit Real, 8-bit offset, Little-endian byte order to load the data
- Type in the width and height of the image
- If you want to open all the files in the folder, select that option



TXM Wizard

- Open TXM Wizard and go to Data Evaluation -> Image Handling
- Click Browse and select the files you would like to load from the file dialog. Make sure the file filter is set to *.binslice
- Click open and view

